

# AI bracelet for health monitoring

Gladkykh Daria, Fateme Hosseini

## **Definition of Heart Sclerosis (Cardiosclerosis)**

- Heart sclerosis is a pathological condition characterized by the replacement of healthy heart muscle (myocardium) with fibrous or scar tissue
- Commonly results from: Myocardial infarction (heart attack), Chronic ischemia, Myocarditis (inflammation), Degenerative aging processes
- This scarring impairs the heart's mechanical and electrical functions, leading to serious clinical manifestations

## **Main Objective**

- To develop an AI-powered system that can detect early signs of heart sclerosis using continuous, non-invasive physiological data collected from a smart bracelet.
- The system will analyze patterns in heart rate variability (HRV), pulse waveform morphology, activity levels, and other biosignals to identify anomalies associated with myocardial fibrosis.

## Secondary Objectives

1. Develop software to collect and stream real-time sensor data from the smart bracelet (T-Watch S3 Plus + MAX30102 PPG sensor)
2. Conduct a pilot study to gather a dataset from volunteers (including healthy individuals and those with known cardiac conditions)
3. Train and validate an AI/ML model capable of flagging potential indicators of heart sclerosis based on wearable-derived metrics
4. Correlate clinical symptoms of cardiosclerosis with measurable biosignals from wearable sensors

## **Relevance & Societal Impact**

- Early detection: Enables timely intervention before heart failure develops
- Accessibility: Non-invasive monitoring for at-risk populations (e.g., post-MI patients, elderly)
- Preventive healthcare: Continuous monitoring supports proactive management of cardiac health

# Symptom–Sensor Correlation Table

Clinical Manifestation of Heart Sclerosis	Relevant Wearable Sensor(s)	Measurable Signal
Arrhythmias / Irregular heartbeat	PPG (MAX30102)	Heart rate variability (HRV), inter-beat intervals
Fatigue, reduced exercise tolerance	Accelerometer + PPG	Activity level, step count, HR recovery post-exertion
Shortness of breath (dyspnea)	Accelerometer + HR trends	Resting HR elevation, abnormal HR response to minimal activity
Conduction abnormalities (e.g., bradycardia)	PPG	Sustained low HR, pauses in rhythm
Fluid retention / nocturnal symptoms	Longitudinal HR & activity trends	Nighttime HR spikes, reduced sleep quality (indirect)

# **What We've Done This Week**

## **Dataset Research & AI Development**

- Conducted extensive research on publicly available PPG, HRV, and Accelerometer datasets.
- Identified a multi-tier strategy prioritizing Foundation Training PulseDB followed by Motion Artifact Handling PPG-DaLiA.
- Began setup and preprocessing of the PulseDB foundation dataset in parallel with hardware integration to maximize model training efficiency.
- Plan to transition to our own collected data once the wrist-worn hardware is operational to fine-tune the final classification model.



# AI/ML DATASET STRATEGY

Category	Dataset	Size / Key Feature	Justification for Use
Foundation Training (Phase 1)	PulseDB (MIMIC-III/VitalDB Derivative)	5,361 subjects / 5.2 Million pre-cleaned 10-s segments of ECG, PPG, ABP.	ML-Ready Baseline: Used to train the initial feature extractor and the <b>first disease classifier</b> using linked MIMIC-III heart failure labels.
Motion Artifact Handling (Phase 2)	PPG-DaLiA + WESAD/WildPPG	136+ subjects with synchronized wrist PPG and 3-axis ACC.	Noise Resilience: Enables the implementation of synthetic augmentation to adapt the baseline model to the high motion noise of our wrist-worn device.
Disease Labels & Validation (Phase 3)	MIMIC-III Derived Cohort + Our Pilot Data	Clinical labels for <b>682 Heart Failure patients</b> and 954 controls.	Ground Truth Proxy: Provides a proven clinical endpoint Heart Failure as the best available proxy for detecting early-stage myocardial fibrosis.

# IMPLEMENTATION ROADMAP: THE 3-PHASE PIPELINE

Phase	Focus & Goal	Key Task & Dataset(s)	Deliverable
Phase 1: Foundation & Classification (Current Focus)	<b>GOAL:</b> Establish the ideal feature extraction and initial disease detection logic.	Train on <b>PulseDB</b> and label a subset using the <b>MIMIC-III Derived HF Cohort</b> .	<b>Initial Disease Classifier &amp; Feature Extractor</b>
Phase 2: Motion Artifact Handling	<b>GOAL:</b> Adapt the Phase 1 model for real-world wrist-worn usage.	Implement <b>Synthetic Augmentation</b> using ACC data from PPG-DaLiA.	<b>Motion Artifact Removal Module</b> (Model is now wrist-robust)
Phase 3: Final Validation & Deployment	<b>GOAL:</b> Test the motion-robust model against custom data and finalize performance metrics.	A) Collect <b>Our Pilot Data</b> . B) Fine-tune and validate against the <b>HF cohort</b> .	<b>Final Binary Classification Algorithm</b> for Heart Sclerosis Risk

## Hardware Progress

### T-Watch S3 Plus - Received & Setup in Progress

- Successfully received the T-Watch devices
- Configuring the built-in accelerometer for activity tracking and movement detection(in progress)
- Integrating external MAX30102 PPG sensor for heart rate monitoring (in progress)
- Sensor placement: Mounted on wrist underside (radial artery position) via I<sup>2</sup>C connection with proper shielding and strain relief for 48-hour continuous wear (only on paper)

## **Power Management Solution**

### **Extended Battery Life Strategy**

- Challenge: Continuous WiFi transmission drains battery too quickly for 48-hour monitoring
- Solution: Local data storage on device memory during the entire monitoring period
- Data transfer: After 48 hours, patient uses WiFi to upload all collected data to server in one session
- Result: Bracelet operates for full 48 hours without recharging, ensuring uninterrupted diagnostic data collection